

Isolating *F. graminearum* from scabby seeds

Potato Dextrose Agar – 1 L

One liter of PDA is enough for approximately thirty 100 x 15mm Petri dishes. Pour 1 liter of water into a large Erlenmeyer flask and add 39 g of PDA. Stir until most of PDA is dissolved. Plug flask with a cotton stopper, cover with aluminum foil, and autoclave at 121°C for 45 minutes, liquid cycle. Let the flask cool slowly to 51°C in an incubator. Pour agar into plates under laminar flow hood. You should wait two days before using the PDA to see if any contaminants have grown on the plates.

Culturing *F. graminearum* isolates

Surface sterilize scabby wheat kernels in 10% bleach solution for 5 minutes. Then, rinse with sterilized water. Place three seeds on each PDA plate. Store plates at room temperature for 2 or 3 days until you see mycelia. Prompt culturing at this point is key to securing clean isolates. Transfer mycelia to a PDA + rifamycin plate (acidified PDA may also be used). Take only one plug from the mycelia tips surrounding each wheat kernel. Each plug is an isolate and should be labeled accordingly. A record should be kept of the source of each isolate. Try to culture 20 – 30 isolates with about 3 isolates from each source. Keep numbers of subcultures to a bare minimum (Dr. Don White, U of I plant pathologist) or isolates may lose their pathogenicity. Seal cultures with parafilm and store in the refrigerator.

Potato Dextrose Agar + Rifamycin – 1 L

As the sterilized PDA is cooling in the incubator, measure out 0.025 g of Rifamycin (antibiotic) and dissolve in 4 ml of EtOH. Aspirate rifamycin solution with a 10 ml syringe add a sterile 0.2 micron syringe filter and dispense solution directly into PDA. Gently swirl agar in large circles as not to create air bubbles.

Note: Rifamycin is light sensitive and must be stored in the dark.

Acidified PDA

Reduce the pH of pre-mixed PDA to around 4.5 using 10% tartaric acid, then autoclave.

Transferring *F. graminearum* onto new PDA (sub-culturing)

- 1) Prepare 2 liters of PDA
- 2) Sterilize hood with 95% EtOH
- 3) Sterilize knife with 95% EtOH and flame
- 4) Cut a 5 mm plug from the hyphal tips of the isolate
- 5) Place plug upside down in the center of the new plate
- 6) Store FHB isolate at room temperature until it has colonized the whole plate
- 7) Check plate for contamination, use for making inoculum or wrap with parafilm, to prevent desiccation, and store in the refrigerator.
- 8) Isolates can remain viable for over a year when stored on PDA in the refrigerator.